

GIALLOY CB / CB-H



Dental

Non – precious dental alloys

CE 0297

Ceramic bonding alloy on Cobalt basis

- **Gialloys** can be melted by torch, centrifuge and vacuum pressure systems
- **Gialloys** are highly corrosion resistant
- **Gialloys** can be lasered

Processing instructions

Modelation

For 100 % flow out of wax copies, minimum wax thickness should be 0.5 mm.

Investing

Gialloys can be casted with most of the investments usual in trade up to max. casting temperature of 1080°C (1980°F).
Recommended investments are **Gilvest HS** for **Gialloy CB** or **Gialloy CB-H**

Casting

Use ceramic crucibles only.
Cast as fast as possible before oxide layer is torn open. Do not force cool or quench the moulds.
Cool down slowly after casting.

Porcelain

Minimum thickness of alloy before firing porcelain is 0.2 mm.
Before applying the ceramic, fire at 960 - 980°C, sandblast with 110 µm aluminum oxide and clean carefully with steamer.
Long term cooling after firing is generally recommended.

Caution

Never put **Gialloys** in acid bath.

All data and recommendations are guidelines only and based on our technical experience. The above recommendations are given to the best of our knowledge. We grant the quality of our products according to our specification. Any further liability cannot be accepted since the proper application of our products is outside our control

GIALLOY CB / CB-H

Typical chemical compositions and physical properties of Gialloy

	Gialloy CB Ceramic bonding alloy on Cobalt basis, Type 5 Nickel and Beryllium free	Gialloy CB-H Ceramic bonding alloy on Cobalt basis Type 4 Nickel and Beryllium free
Norm	EN ISO 22674	EN ISO 22674
Co	61	~ 58
C	< 0,1	0,04
Cr	28	25
Si	1,7	1
Mo	-	1
W	8,5	9
Fe	< 0,5	1
Mn	0,25	-
Nb		2
V		2
Hardness: (HV 10)	285	340
Elongation: (%)	10	5
Tensile strength: (MPa)	845	790
0,2% Elastic limit: (MPa)	620	
Modulus of elasticity (GPa)	ca. 190	
Coeff. of thermal Expansion (20° - 600° C) (10 ⁻⁶ /K)	14,1	14,1
Melting range (°C) (Solidus/Liquidus)	1390 – 1415	1370 – 1385
Pre-heating temp. of ring (°C)	910 – 940	900 - 930
Density (g/cm ³)	8,3	8,4
	Recommendation: Before applying the ceramic, fire at 960°C	Recommendation: Before applying the ceramic, fire at 980°C