

Gilvest Liquid "G"



Dental

Allows increased expansion values of investments for special requirements, e.g. non precious alloys and telescoping works.

For detailed working instructions please pay attention to the data sheets of the various GILVEST products. This Gilvest liquid **G** raises the value of the stated setting expansion by 35 - 40 %.

GILVEST LIQUID **G** concentrated freezes at temperatures below 0°C. Once frozen it is unusable. Therefore fill your winter requirements early.

Setting expansion of Gilvest Liquid and Gilvest Liquid **G**

